

AI-Enhanced Smartphone-Based Fluorescence Microscopy with IoT-Enabled Staining for Point-of-Care Diagnostics




Transparent plexiglass BLIPS Ultra lens Slide with liver tissue sections Light source





What we Know so far







The BLIPS lens-based imaging system, a The BLIPS lens-based unphone microscopic system along with the acquired images of humiliver tissue section obtained and comparative image of the bright-d microscope. Reproduced from Cesaretti et al. [20], with permission

from Jhon Wiley and Sons (copyright 2016). **b** The real-tracquisition of cancer cells. Reproduced from Skandarajah et al permission from PLOS (copyright 2017)

Working of compound microscope






Types of Mobile microscopes



MICROSCOPE >	AMSCOPE SM-4TZ-144A	AMSCOPE ME580TA-PZ-2L	SHIIRI	OXBIRD	ANNLOV	CAINDA	BYSAMEYEE 4KHD WIFI	TOMLOV DM201
ADVERTISED MAGNIFICATION	90-225X	500X	1000X	2000X	1000X	1000X	1000X	1200X
ACTUAL MAGNIFICATION	225X	500X	~100X	~100X	~100X	~100X	~50X	~50X
SENSOR	N/A	N/A	инкножн	инкножн	GC0308	GC0308	UNKNOWN	UNKNOWN
ADVERTISED RESOLUTION	N/A	N/A	NOT SPEC'D	NOT SPEC'D	8MP 1280X720	NOT SPEC'D	8MP 3840X2160	16MP
ACTUAL RESOLUTION	N/A	N/A	O.3MP 640X480	O.3MP 640X480	0.3MP 640X480	0.3MP 640X480	8MP 3840X2160	16MP 3840X2160
CONNECTIVITY	N/A	N/A	USB	USB	USB	USB	WIFI & USB	MICRO SD, HDMI
PRICE	\$654.99	\$992.99	\$9.99	\$12.99	\$49.99	\$22.99	\$39.99	\$169.99
SONIC DOODLE (NORMALIZED)	ES V							Elli.

Types of Digital microscopes



Types of Compound microscopes






Types of Microscope cameras

Limitations

- Lighting
- Design



Development of a low-cost smartphone fluorescence microscope.







Figure 2. Use of recreational LED flashlights and theater stage lighting filters for smartphone green fluorescence microscopy. (a) Schematic of components used for green fluorescence viewing on the glowscope.


Fig. 4 Smartphone-based fluorescence microscope. A The optical design of the smartphone fluorescence microscope. B System photo and design sketch of the smartphone fluorescence microscope. (Figure reproduced from Shan et al., 2019 with permission from Elsevier).





158 Biophysical Reviews (2024) 16:151–163

Fig. 5 a 3-D model of the threewavelength illumination unit and b the mobile prototype with a smartphone on it. (1) Laser modules (three pairs, 448–532– 659 nm). (2) shielding cylinder, (3) beams collector, (4) diffuser of LASER light, (5) sticky platform, and (6) electronics compartment. (Figure reproduced from Spigulis et al., 2017 with permission from SPIE).



a. Traditional Method of Bacteria Identification



b. BIOINTEL Proposed Method of Bacteria Identification




Fig. 1. Comprehensive pipeline of BioIntel system.

Cell/Tissue analysis





Fig. 6 Applications to imaging cells with various illumination patterns. a Diatoms and E. gracilis images zoomed in and cropped under bright-field (BF), dark-field (DF), and Rheinberg (RI) illumination patterns (second and third rows). In the top row, these illumination patterns were created using the Retina display. Used ring illumination patterns with single dark-field (DF-single) and double dark-field (DFdouble) (second and third columns). On a dark specimen background, different types of RI were prepared: blue-yellow patterns (fourth col-


umn), and half-rings with multiple colors (red-green half-ring pattern) (fifth column). **b** Images of stained (methyl blue) and unstained human epithelial cheek cells (HECC) under BF and DF illumination patterns, zoomed in and cropped. Images from the second column were subtracted using a single dark-field (DF-single) illumination pattern. (Figure has been reproduced from Kheireddine S. et al., 2019 with permission from RSC publication).



g. 3 Histology images acquired with Pocket MUSE. All samples were stained with 0.05% w/v Rhodamine B and 0.01% w/v DAPI unless otherwise ecified. a Image of a thick section of mouse kidney sliced with a razor blade. A close-up view of the region in the white box (box size: 320 × 320 µm²) is own on the right. b Image of mouse skeletal muscle torn with tweezers. c Image of the serosal surface of a mouse small intestine. d Image of salmon eak sliced with a razor blade. An additional 0.1% w/v Light Green SF dye was added

Bacterial Sample Identification

Limitation of resolution






Fig. 9 Images of bacterial samples acquired with Pocket MUSE. Two different populations (Gram + and Gram +) of bacteria can be separated with Pocket MUSE. Aliquots of (a) deionized water, (b) an *E. coli* (Gram +) culture, (c) a *B. subtilis* (Gram +) culture, and (d) a mixture of *E. coli* culture and *B. subtilis* culture imaged with Pocket MUSE. All samples were stained with 0.01% w/v DAPI and 0.02 % w/v WGA-AF594, and imaged under the same conditions. Images show (a) no bright spots, (b) blue-grayish spots, (c) a mixture of reddish and orangish spots, and (d) a mixture of blue-grayish, reddish and orangish spots dispersed across the FOV. For each image, a close-up view of the region labeled with the box (box size: $250 \times 250 \,\mu\text{m}^2$) is shown below. For (b)-(d), the center pixels of three bright spots in each image, pointed to by the arrows, are shown in the boxes below. The RGB values (top to bottom) of the pixels are labeled on the side. The nucleic acid in *E. coli* (DAPI stain) contributes to the blue-grayish speckles ($R \approx G \approx B$). Peptidoglycan on the *B. subtilis* surface (WGA-AF594 stain) contributes to the reddish speckles (R > G + B). Nucleic-acid-rich endospores in *B. subtilis* contribute to the orangish

Proposal

Smartphone + AI + IoT = Point-of-Care Diagnosis

Smartphone fluorescence microscope

QR-coded slides → automatic model + staining selection

IoT-enabled staining system (CISS + Raspberry Pi)

Al-powered detection and diagnosis

Key Modules

- Smartphone Microscope with fluorescence
- Al for:
- Image preprocessing
- Pathogen classification
- Quality check
- IoT-Controlled Staining Unit:
- Raspberry Pi/Arduino + CISS
- QR-code-triggered stain protocol



Applications

- Infectious disease detection: malaria, TB, typhoid in rural area
- Water quality testing (pathogen load)
- Blood/sputum diagnostics in remote clinics
- Telemedicine integration
- Food Quality

School/university science labs

Innovative Approach to Diagnosis

- Smartphone + AI + IoT = Point-of-Care Diagnosis
- Smartphone fluorescence microscope for rapid analysis
- QR-coded slides enable automatic model and stain selection
- IoT-enabled staining system using CISS and Raspberry Pi
- Al-powered detection ensures accurate diagnosis